
	ITEM NUMBER:	DRAWN BY: BVD	TITLE: Chip Serpentine 140x250um, 500mm			
	THIRD ANGLE PROJECTION 	DRAWN DATE: 04 Jul 2025				
		PROJECT NO: RD0946.001				
<p><u>PROPRIETARY AND CONFIDENTIAL</u></p> <p>THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF MICRONIT BV. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF MICRONIT BV IS PROHIBITED.</p>	UNLESS OTHERWISE SPECIFIED:	EXT. REF.:	DWG NO. 42004744		REVISION 2.0	
		MATERIAL: Glass D263 Eco				
	- DIMENSIONS ARE IN MILLIMETERS - LINEAR TOLERANCES: ± 0.1 - ANGULAR TOLERANCES: ± 5°	SCALE: 2:1		SHEET FORMAT: A4		SHEET 1 OF 1
	DO NOT SCALE DRAWING					