



Dimension throat	
After etching	50 μm
Mask width	10 μm
Depth	20 μm
Dimension pore	
After etching diameter	90 μm
Mask width diameter	50 μm
Depth	20 μm

Note:
The uniform network consists of 10151 pores, which are connected by throats.
In this number is included the pores connected to the inlet/outlet channels

Zone	Rev	REVISION NOTE	Mod By	App By	App Date
	1.1		IR		05 Jun 2023
	1.1	Conformity changes	ME		
	1.0	First Release	ME		09 Aug 2022



ITEM NUMBER: 02978
DRAWN BY: WVDC
DRAWN DATE: 22 Dec 2021
PROJECT NO:
EXT. REF.:

TITLE: Chip EOR UN 20.2

PROPRIETARY AND CONFIDENTIAL	UNLESS OTHERWISE SPECIFIED: - DIMENSIONS ARE IN MILLIMETERS - LINEAR TOLERANCES: ± 0.1 - ANGULAR TOLERANCES: ± 5°	MATERIAL: Borofloat 33	DWG NO. 42005552	REVISION 1.1
THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF MICRONIT BV. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF MICRONIT BV IS PROHIBITED.		DO NOT SCALE DRAWING	SCALE: 5:1	SHEET FORMAT: A3 SHEET 1 OF 1