



DETAIL B
SCALE 50 : 1



SECTION A-A
SCALE 5 : 1



SECTION C-C
SCALE 50 : 1

The structures are randomly placed, resulting in throats and channels of random dimensions. The design resembles structures from cut sandstone.

Zone	1.0	First Release	Mod By	IR	05 Jun 2023
Rev		REVISION NOTE	App By		App Date

	ITEM NUMBER:	02976	DRAWN BY:	ME	TITLE:	Chip EOR PR 20.2		
	THIRD ANGLE PROJECTION		DRAWN DATE:	11 Nov 2022	DWG NO.	42005587	REVISION	1.0
PROPRIETARY AND CONFIDENTIAL	UNLESS OTHERWISE SPECIFIED:	- DIMENSIONS ARE IN MILLIMETERS - LINEAR TOLERANCES: ± 0.1 - ANGULAR TOLERANCES: ± 5°	MATERIAL:	Borofloat 33	SCALE:	5:1	SHEET FORMAT:	A3
THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF MICRONIT BV. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF MICRONIT BV IS PROHIBITED.	DO NOT SCALE DRAWING							SHEET 1 OF 1