



SHT-G.025-01



2 x  $\phi 0.40 \pm 0.10$

SECTION B-B

SECTION C-C

SECTION D-D

DETAIL E  
SCALE 50 : 1

SECTION F-F  
SCALE 50 : 1

DETAIL A  
SCALE 50 : 1

	1.0	First Release		NVDB	13 Oct 2023
Zone	Rev	REVISION NOTE	Mod By	App By	App Date



ITEM NUMBER: 02134  
 DRAWN BY: BVD  
 DRAWN DATE: 13 Oct 2023  
 PROJECT NO: RD0945.001  
 EXT. REF.:

TITLE: Chip glass 2-port sheet nozzle

PROPRIETARY AND CONFIDENTIAL  
 THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF MICRONIT BV. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF MICRONIT BV IS PROHIBITED.

UNLESS OTHERWISE SPECIFIED:  
 - DIMENSIONS ARE IN MILLIMETERS  
 - LINEAR TOLERANCES:  $\pm 0.1$   
 - ANGULAR TOLERANCES:  $\pm 5^\circ$   
 DO NOT SCALE DRAWING

MATERIAL: D263  
 DWG NO. 42008759  
 REVISION 1.0  
 SCALE: 10:1  
 SHEET FORMAT: A3  
 SHEET 1 OF 1